

# Image analysis software

<https://github.com/albagranados/cellviewer/tree/master/imageanalysis>

Alba Granados

Jérôme Solon's laboratory  
Center for Genomic Regulation, Barcelona, Spain

October 10, 2018



# Outline

1 Image pre-processing

2 Image processing

3 Image analysis

## Workflow 0: read dataset

```
file_dirs = [ '/home/alba/ownCloud/postdoc_CRG/coding/github/
.....cellviewer/data/melike/Histones/crops/noTSA/',
              '/home/alba/ownCloud/postdoc_CRG/coding/github/
.....cellviewer/data/melike/Histones/crops/TSA/']
is_dataset = 0
file_name = 'hFb.2_driftcorrectedlist_85120_105135 '
fileExt = '.txt'
is_storm = 0
run = dict(image_processing=1, image_analysis=1)
```

**file\_dirs** contains one or more paths to the experiments

**is\_dataset** is 0 if we only want to analyze 1 cell/window

**file\_name**, fill if is\_dataset=0,

**is\_storm** is 1 if the file contains all output columns from typical STORM output. It is 0 if the .txt file is two-columns with x-corrected and y-corrected (e.g., crop window from entire STORM image)

## Workflow 0: bin2txt and crop

Input parameters:

```
dict_inputfile = dict(file_dir=file_dir , file_name=file_name ,
                      fileExt=fileExt , is_storm=is_storm ,
                      out_channel='all' ,
                      ispp=1, compute_ROI=0, crop=1,
                      crop_range=[85, 120, 105, 135] ,
                      pixelate=0,
                      tessellate=1,
                      original_pixel_size=160, photonconv=0.14,
                      resolution=1)
```

**out\_channel** matters if **is\_storm**=1. Values: 'all', [1,2], [0],....

**ispp** is 1 if input is a point pattern, like STORM output, and not a regular image.

**compute\_ROI** is 1 if we want to compute a smaller area. Then, if **crop** is 1, we can select the **crop\_range** as  $[x_0, x_1, y_0, y_1]$ .

**pixelate** is 1 if the image is generated by regular pixellation, or **tessellate** is 1 if we do that via Voronoï tessellation.

**original\_pixel\_size** is STORM pixel size (nm)

# Outline

1 Image pre-processing

2 Image processing

3 Image analysis

## Workflow 1: Point pattern to image

Input parameters:

```
analysis_pixel_size = 10
scale_pixel_size = float(analysis_pixel_size)/
                    dict_inputfile.get('original_pixel_size')
dict_image = dict(scale_pixel_size=scale_pixel_size,
                  original_pixel_size=dict_inputfile.
                    get('original_pixel_size'),
                  interpolate_method='linear',
                  detect_densitytransform='log',
                  descr_densitytransform='linear')
```

**analysis\_pixel\_size** is the pixel size (regular grid) of the analysis density-based image (point(STORM) to image).

**interpolate\_method** is the method of density interpolation from irregular grid (Voronoi) to regular grid.

The Voronoi-based density map (image) is transformed for the feature detection algorithm (**detect\_densitytransform**) and the feature descriptor algorithm (**descr\_densitytransform**); go to slide 8

## Workflow 2: Cluster detection

Input parameters:

```
print '\n-----IMAGE_PROCESSING-----'
dict_sift = dict(scale_pixel_size=scale_pixel_size, resolution=1,
                 original_pixel_size=dict_inputfile.
                                     get('original_pixel_size'),

                 # feature detection
                 t=10, feature_name='blob', # 0.8, scale_ini=80
                 thresholding=1, threshold_percent=0.8,
                 num_features='all', scale_range_is='nm', scale_ini=40,
                 scale_end=200,
                 scale_spacing='odd', nscales=150,
                 scale_resolution=dict_inputfile.get('resolution'),
                 max_filter_width=7, max_filter_depth=7,

                 [...])
```

**t** is a single scale, but typically we'll use a range (see below)

**thresholding** can be activated. **threshold\_percent** corresponds to the % of the largest values of the Laplacian.

**num\_features** is the maximum number of features to be detected. Can be 'all'

**scale\_ini** and **scale\_end** defines the limits of the scale range search. If 'odd' is

**scale\_spacing**, then  $t_1, \dots, t_n$  s.t.  $3\sigma = 3\sqrt{t_i} = d_i$  for  $d_i = 2n + 1$  diameters (odd number of pixels), and **nscales** is neglected.

Filter for local maxima has dimensions **max\_filter\_width** (space) and **max\_filter\_depth** (scale)

## Workflow 2: Cluster detection

### `detect_densitytransform`





Figure: Detection with a logarithmic (left) and linear (right) density map. Logarithmic is shown in both cases for clarity reasons.

### `threshold_percent` needed?



Noise is discernible from circular clusters in terms of the Laplacian.



## Workflow 3: Cluster description

SIFT input parameters:

```
dict_sift = dict(  
    [...]  
  
    # feature description [main orientation(s)]  
    compute_orientation=1, n_bins_ori=36, peak_ratio=0.7,  
    sigma_ori_times=1.5, window_ori_radtimes=1,  
    smooth_cycles=2,  
  
    # feature description [histograms]  
    compute_sift_descr=1,  
    sigma_descr_times=2, window_descr_radtimes=1,  
    n_hist=16, n_bins_descr=8, threshold_sat=0.2,  
    plot_graphics=0)
```

**n\_bins\_ori** is the number of bins of the histogram of gradients arguments to define the main orientations.

**peak\_ratio** is the portion of the largest bin (largest orientation) above which all maxima are considered main orientations.

**smooth\_cycles** is the number of cycles to smooth the histogram.

**sigma\_ori\_times** and **window\_ori\_radtimes** define the window of the local patch to compute the gradients. Recall  $\sigma = \sqrt{t}$  (scale).

**sigma\_descr\_times** and **window\_descr\_radtimes** is the size of the local patch to compute the  $n$ -dimensional SIFT descriptor, where  $n = \mathbf{n\_hist} \cdot \mathbf{n\_bins\_descr}$ .

With **threshold\_sat** we reduce the influence of large gradient magnitudes by thresholding the values in the unit feature vector to each be no larger than 0.2, and then renormalizing to unit length.

# Outline

1 Image pre-processing

2 Image processing

3 Image analysis

## Workflow 4: Cluster-based classification

```
print '\n-----IMAGE_ANALYSIS-----'; ini_time = time.time()
init = 'k-means++'
k0 = 2; kn = 5; serror = [] # build bag of words
for k in range(k0, kn+1):
    [...]
```

**k0** and **kn** define the range of number of clusters/words we want to analyze